

Features

- Class AAA specification (ASTM, IEC)
- Illumination area: 50x50mm
- Touchscreen power supply with control software included
- Manual shutter included (electronic shutter available)
- Variable attenuator from 0.1–2 suns
- Plug and play operation
- Long working distance can facilitate glovebox integration

Applications

- Photovoltaic Testing
- Environmental Testing
- Photobiology and Photochemistry
- Material and degradation testing





Overview

Sciencetech's line of SciSun solar simulators are easy to use, economically priced, and technically superior. The SciSun line is designed for researchers who do not require a large field of illumination. They can produce up to **2 Suns** and feature Class AAA specifications.

The SciSun series provides a flexible output orientation that can be adapted to different requirements. The standard configuration is downward-facing; however, a horizontal output can be achieved easily.

All SciSun models include:

- arc lamp housing with integrated igniter
- xenon arc lamp
- filter holder
- beam turner (variable illumination directions)
- quality control report

Non-LP series models also include:

- touchscreen power supply interface
- power supply control software
- manual variable attenuator
- height-adjustable stand

STANDARDS

SciSun solar simulator specifications listed are according to ASTM E927-19 and IEC-60904-9 unless otherwise stated.

Specifications

Model	SciSun-300	SciSun-LP-300	SciSun-150	SciSun-LP-150
Part Number	160-9101	160-9104	160-9103	160-9105
Target Area	50 × 50 mm			
Irradiance Uniformity	Class A ¹			
Irradiance at Target (AMI.5G Sun=100mW/cm²)	Up to 2 Sun ²		Up to 1 Sun ²	
Lamp Wattage (watts)	300		150	
Spectral Match AM 1.5G	Class A ³			
Lamp Type	Xenon Short Arc , Ozone free			
Temporal Stability	Class A ⁴			
Working Distance (mm)	380 ± 15			
Manual Shutter	Included			
Manual Variable Attenuator	Included	Available	Included	Available
Dimensions (L×W×H)	535 × 183 × 188 mm			
Weight without PS (kg)	8.5 + 8 (stand)	8.5	8.5 + 8 (stand)	8.5
Power Supply Model	601-300	EPS-300	601-150	EPS-150
Power Requirements	110-240V, 50Hz/60Hz , 450W			
Stability / Ripple / Regulation	0.05% / < 1% / 0.02% current variation for 5V line charge			

¹⁾ Determined from true Isc measurements with silicon sensor mounted on 2 axis automated stage. 2) Measured using NIST traceable secondary reference cell

³⁾ Measured with scanning spectroradiometer calibrated as per ASTM G138-06 4) Determined from 20 measurements spaced at 250ms, NPLC=1







Solar Simulator Classification Measurement

Class A Spectral Match Measurement:

SciSun solar simulators match Class A spectral match when used with a compatible air mass filter (sold separately; see below using an AMI.5G filter). All testing results are for an example SciSun-150 and individual reports will vary.






Class A spatial non-uniformity (NU):

SciSun solar simulators meet Class A spatial non-uniformity by default (see below).

Non-uniformity = **1.3**%


Class B may also be available over larger target sizes upon request.



Class A Temporal Instability:

SciSun solar simulators meet Class A temporal instability.

Temporal Instability of SciSun-150



STANDARDS

SciSun solar simulator specifications listed are according to ASTM E927-19 and IEC-60904-9 unless otherwise stated. We can accommodate testing to match several standards.

Testing procedure as per ASTM E927-19 provided by default. Please specify upon ordering if testing against IEC-60904-9 is required.



Standard Features

FILTER BOX ASSEMBLY

Spectral Filter Options			
Model	Description		
AMI.5G-FT-3	AMI.5G Filter—Class A		
AMI.0D-FT-3	AMI.0D Filter—Class A		
AMI.5D-FT-3	AMI.5D Filter—Class A		

This system has a modular optics assembly which can hold a range of filters in Sciencetech's standard FT style filter holder. The most popular options are AM filters; however a range of other filter options are available such as bandpass filters and neutral density filters.

Variable Aperture VAR-ATTN-M

Sciencetech's SciSun solar simulators include a variable aperture component, which allows variation of the output irradiance level without adjusting the power supply. The range of attenuation is continuously variable from 10% to 100%. Uniformity is best maintained at specific output levels.

Non-uniformity versus output level for the VAR-ATTN-M may vary between models.

POWER SUPPLY AND SOFWARE CONTROL

Each SciSun series solar simulator (non-LP series) comes with a 601-series power supply.

Standard features included with Sciencetech's 601–series power supplies:

- Touchscreen interface
- Shutter and exposure control*
- Single connection for lamp power, cooling, and communication
- Lamp starts and timer log
- Fan cooling safety interlock
- RS232 software GUI included

Control 1.1

Control Shutter © Open © Close
Cooling © ON © OFF
Lamp OFF
Lamp OFF
Current 0.0

Voltage 0.0

Voltage 0.0

Power 0

Refresh Feedback Status

Refresh Comm

Microsoch NET 4.5.2 Framework programming by Jim Sample

COMS

SCIENCIFICAL

Shutter Closed
Cooling OFF
Lamp OFF
Current 0.0

Voltage 0.0

Power 0

Total Lamp Starts 0

Total Lamp Timer 0 Hours 0 Min

Software GUI for power supply control

The SciSun-LP series features the EPS-series simplified power supply, which lacks a touchscreen, computer control, or control of electronic accessories.




601-series power supply

601-series touchscreen power supply main control screen



601-series touchscreen power supply automatic shutter control screen







Popular Accessories



Dimensions

Dimensions are in [mm] and inches. [406.4] 16.00 [615.9] 24.25 [711.2] 28.00 [670.6] 26.40 6.4 SciSun-150 and SciSun-300 SciSun-150 and SciSun-300 675 ROTATION 360° [188] 7.38 SciSun-LP-150 and SciSun-LP-300 [183] 7.20

